


BARIATRIC TOILET SEAT RAISER



**MAX USER WEIGHT
160kg**

TOILET NOT INCLUDED

**Width between handles
500mm**



FEATURES

RG520

- Assist user when toileting by elevating the toilet seat height
- Armrests to help getting on and off the toilet
- Height adjustable seat
- Plastic flip up toilet seat flap and lid
- Removable splashguard
- Non-slip leg tips
- Flat pack design
- Tool free assembly
- Can be converted to a commode
RG5B: Bowl and lid (***sold seperately***)

Please turn over for product specifications



SCAN THE QR CODE TO
CHECK OUT THE FULL
REDGUMBRAND RANGE!



**COMPLIANCE
TESTED**
AS/NZS3973-2009


ASK FOR THE BRAND

redgumbrandmobility.com.au


SPECIFICATIONS

(All specifications may vary (+/- 10mm))


Overall measurements are taken from the lowest seat height setting




Overall height
650mm (H)




Width between armrests
500mm (W)



Seat width
400mm (W)



Seat opening
180mm (W)
230mm (W)
275mm (L)



Seat height
(floor to top of seat
flap when closed)


Hole 5: **545mm**

Hole 4: **520mm**


Hole 3: **495mm**

Hole 2: **470mm**


Hole 1: **445mm**




Overall width
625mm (W)




Minimum toilet clearance
385mm (H)




Seat depth
400mm (D)



Splash guard depth
140mm (D)



Overall depth
495mm (D)



Armrest height
200mm (H)



Product weight
5kg



Maximum user weight
160kg


Benefits can be achieved for some users by
adjusting the angle of the toilet seat raiser.

Seat tilt settings

2 hole variation = **96°**

1 hole variation = **93°**

No variation = **90°**



WARNING

Front legs must
stay at 90° to the
floor at all times.

Do NOT have
more than a 2 hole
variation between
the front and rear
legs. The maximum
angle achievable
is 96°.

ASK FOR THE BRAND

redgumbrandmobility.com.au